



## Product Datasheet

### Allergy Panel – Specific IgE

### True Human Serum off the Clot

#### Application:

The in.vent Panels of Allergy are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. They are representative for different Allergens. The entire Range of Laboratory Test Results is represented within the Panels. The Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

<b>Product No.:</b>	ALPA 4.1.7.3.1
<b>Number of Samples</b>	10
<b>Volume per Samples</b>	1.0 ml
<b>Source Material</b>	true human serum off the clot
<b>Details given</b>	product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life
<b>Biomarker for Choice</b>	CCD, d1, d2, e1, e2, e3, e204, f1, f2, f3, f4, f5, f9, f10, f13, f14, f17, f20, f23, f24, f25, f27, f31, f35, f45, f49, f75, f76, f77, f78, f84, f85, f88, f237, g1, g3, g6, g12, h1, i1, i3, i6, k82, m1, m2, m3, m6, t2, t3, t4, t7, t23, w1, w6, w9
<b>Provided Titer</b>	2x normal, 4x low positive, 4x high positive
<b>Disease State Testing</b>	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV, RNA, Syphilis.
<b>Materials Safety</b>	These methods cannot eliminate possible infectivity of the material. All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material.  All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	24 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)