



## Product Datasheet

### Pool of Human Citrate-Phosphate-Dextrose (CPD) Plasma

#### Recovered Human Plasma (from whole blood)

#### Application:

Human CPD Plasma is designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. Human CPD Plasma is used for production purposes, like production of standards, calibrators or controls for IVDs or EQAS. For non-therapeutic, in-vitro-diagnostic use only.

<b>Product No.:</b>	NDPO 4.2.1.2
<b>Bottling</b>	100 ml
<b>Source Material</b>	whole blood plasma (CPD) from healthy donors, frozen within < 96 h, non-lipaemic, non-haemolytic
<b>Filtration</b>	0.2 µm
<b>Details given</b>	product number, volume, matrix, minimum shelf life, result of Laboratory testing (if applicable)
<b>Test Samples</b>	1.0 ml
<b>Disease State Testing</b>	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV, RNA, Syphilis.
<b>Materials Safety</b>	These methods cannot eliminate possible infectivity of the material. All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material.  All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	24 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](http://inquiry.inventdiagnostica.de)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)