

Product Datasheet

Human Thyroid Tissue

Application:

Human Tissues are a high potential resource for extraction of many different proteins. Fresh frozen Tissues are essential for biomarker research and successful in-vitro diagnostics development. All Thyroid Tissues provided by in.vent Diagnostica are suitable for extraction of various kinds of Human Thyroidperoxidase (TPO).

Product No.:	HTT 5.1.1
Volume	> 5 gram per donor, individual weight on request
Source Material	surgically removed Thyroid Tissue because of medical indications, tissue frozen at -20°C within 2 hours after resection, region of origin: Central and Eastern Europe
Titer	not applicable
Details given on Certificate	weight, viral testing, minimum shelf life
Informed Donor Consent	existing
IRB/ IEC	existing
Test Samples	not available
Viral Testing	Each individual tissue was partly tested for HIV-1/2 RNA, HCV RNA and HBV DNA with CE marked test systems and found negative or non-reactive. These methods cannot guarantee the absence of an infectious agent. Therefore standard precautions as if handling potentially infectious material are advised.
Materials Safety	All products are for research or manufacturing purposes only and are not for use in human therapeutic procedures.
Shipping	on dry ice
Storage	≤ 18°C
Shelf Life	48 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Central Contact

+49 (0)3302 55 199 0
info@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Contact for your Order

sales@inventdiagnostica.de

