



Product Datasheet

Lyme Disease Panel

True Human Serum off the Clot or Plasma

Application:

The in.vent Panels of Lyme Diseases are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. They are representative for different types of Borrelia Species. The entire Range of Laboratory Test Results is represented within the Panels. Laboratory reports include LINE Blot and/or ELISA results. The Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

| | |
|-------------------------------------|--|
| Product No.: | DSPA 4.1.4.3.1 |
| Number of Samples | 10 |
| Volume per Samples | 1.0 ml |
| Source Material | true human serum off the clot or plasma |
| Details given | product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life |
| Disease State Testing | Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV, RNA, Syphilis. |
| Materials Safety | These methods cannot eliminate possible infectivity of the material. All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals. |
| Storage | ≤ -15°C |
| Minimum Shelf Life | 24 months |
| Manufacturing Quality System | in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485. |
| Pricing | Please inquire for your individual quote. |

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de