



Product Datasheet

In.QualiSet Sepsis Panel

Sepsis associated serum and plasma

Application:

The diagnosis of sepsis remains a major challenge, which is to be overcome with our range of samples. The In.QualiSet Sepsis Panel is specially designed for in-vitro diagnostic test kits for research through to quality control, targeted on the promising sepsis biomarker PCT.

Product No.:	DSPA 4.1.5.7.1
Number of Samples	96
Volume	150 µL per sample
Matrix	serum and EDTA plasma
Product Details	product number, sample number, volume, matrix, age, gender, measured laboratory result, patient documentation, medication, diagnosis, assay system, minimum shelf life
Product Content	<p>samples with sepsis diagnosis:</p> <ul style="list-style-type: none"> • 8 x progress samples (measured on PCT, Bilirubin, CRP, Creatinin, D-Dimer, Lactat) • 32 x PCT low positive ($\geq 0,5 < 2$ ng/ml) • 32 x PCT positive ($\geq 2 < 10$ ng/ml) • 8 x PCT high positive (≥ 10 ng/ml) • 8 x serum-plasma duo PCT higher positive ($\geq 0,5 < 10$ ng/ml)
Infectiology	Samples of human origin are potentially infectious. Therefore, the necessary safety measures should be taken when handling it.
Application	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	$\leq -15^{\circ}\text{C}$
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual offer

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de