



Product Datasheet

Toxoplasma gondii IgM

Human Plasma Pool, anti-Toxoplasma gondii IgM positive

Application:

Pool of Human Plasma, anti-Toxoplasma IgM positive, is used for production purposes like production of standards, calibrators or controls for IVDs or EQAS.

Product No.:	DSPO 4.1.7.1.1
Volume	10.0 l individual volumes on request
Source Material	true human CPD plasma, pooled
Titer	2.10 COI
Details given	titer, volume, matrix, test kit, test manufacturer
Bottling	1.0 l
Test Samples	available
Disease State Testing	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV-I NAT, HCV NAT.
Materials Safety	All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	48 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de