



## Product Datasheet

### Calprotectin-Elastase Control

#### Faeces

#### Application:

in.vent Quality Controls are designed to be used in Laboratory Quality Control Settings, in Quality Departments and R&D Units. They are applicable for solving tasks of precision, accuracy or correctness of in-vitro diagnostic assay measurements and analyte components throughout the whole value chain of in-vitro diagnostic test kits or assays.

Faecal Calprotectin is used to detect intestinal inflammation, especially Ulcerative Colitis, Crohn's Disease and other Inflammatory Bowel Diseases.

<b>Product No.:</b>	QC 2.1.12
<b>Source Material</b>	faeces
<b>Bottling</b>	0.5 g each vial
<b>Provided Levels</b>	calprotectin healthy, elastase-pathogen: <ul style="list-style-type: none"> <li>level: 0 µg/g calprotectin, ~15 µg/g elastase</li> </ul> calprotectin-pathogen, elastase healthy: <ul style="list-style-type: none"> <li>level: ~900 µg/g calprotectin, ~400 µg/g elastase</li> </ul>
<b>Stability</b>	12 months at 2-8 °C
<b>Form</b>	lyophilised
<b>Test System</b>	calprotectin: fCALTM ELISA Calprotectin, Bühlmann elastase: Pancreatic Elastase 1 ELISA, ScheBo Biotech AG
<b>Disease State Testing</b>	on request (TPE-test for salmonella, yersinia, shigella and campylobacter)
<b>Pre-Samples</b>	available
<b>Delivered Information on Certificate</b>	product number, matrix, processing, homogeneity, target values, reconstitution instruction
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)