



## Product Datasheet

### Diabetes Type I & II Panel

#### True Human Serum off the Clot

#### Application:

The in.vent Panels of Diabetes are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. The entire Range of Laboratory Test Results is represented within the Panels. The Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

<b>Product No.:</b>	DSPA 4.1.3.3.1
<b>Number of Samples</b>	20
<b>Volume per Sample</b>	1.0 ml
<b>Source Material</b>	true human serum off the clot
<b>Details given</b>	product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life
<b>Diagnoses</b>	<ul style="list-style-type: none"><li>• 10 samples of Diabetes Type I (measured values on availability: anti-GAD, IA-2 antibodies, IAA, anti-Zinc transporter 8)</li><li>• 10 samples of Diabetes Type II (measured values on availability: HbA1c, Blood glucose)</li></ul>
<b>Materials Safety</b>	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	24 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](http://inquiry.inventdiagnostica.de)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)