



Product Datasheet

Positive SARS-CoV-2 Panel

True Human Serum off the Clot

Application:

The in.vent Positive SARS-CoV-2 Panels are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

The entire Range of Laboratory Test Results is represented within the Samples. The

Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

| | |
|-------------------------------------|--|
| Product No.: | DSPA 4.1.9.3.1 |
| Number of Samples | 8 |
| Volume per Samples | 1.0 ml |
| Source Material | true human serum off the clot |
| Details given | product number, sample number, matrix, volume, age, gender, symptoms start date, date of positive PCR, titer (IgG and IgM), test system, manufacturer of test system, minimum shelf life, date of blood donation |
| Diagnoses | 2 serum IgG low positive/IgM negative (= 2 samples) 1 serum sample IgG clear positive/IgM clear positive (= 1 sample) 1 serum sample IgG clear positive/IgM negative (= 1 sample) 2 consecutive samples (two donations IgG and IgM positive) (= 2 samples) 2 consecutive samples (both clear positive and significantly decreasing for IgG and IgM), (= 2 samples) |
| Materials Safety | All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals. |
| Storage | ≤ -18°C |
| Minimum Shelf Life | 24 months |
| Manufacturing Quality System | in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485. |
| Pricing | Please inquire for your individual quote. |

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de