



Product Datasheet

Prostate Specific Antigen (PSA & fPSA) Control, CE marked

True Human Serum off the Clot

Application:

in.vent Quality Controls are designed to be used in Laboratory Quality Control Settings, in Quality Departments and R&D Units. They are applicable for solving tasks of precision, accuracy or correctness of in-vitro diagnostic assay measurements and analyte components throughout the whole value chain of in-vitro diagnostic test kits or assays.

Product No.:	QC 2.1.10
Source Material	true human serum off the clot
Bottling	0.5 ml each vial
Provided Levels	level 1: 3-6 ng/ml PSA, 1-3 ng/ml fPSA level 2: 6-9 ng/ml PSA, 3-6 ng/ml fPSA
Stability	1 year at 2-8 °C
Form	lyophilised
Disease State Testing	tested by CE marked test systems and found negative for HIV I/II antibodies, HBs antigen, HCV antibodies
Pre-Samples	available
Delivered Information on Certificate	product number, matrix, target values, assay system, minimum shelf life, homogeneity, result of infectious disease testing, reconstitution instruction
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de