



Product Datasheet

Normal Pregnancy Panel – Consecutive Donations

True Human Serum off the Clot

Application:

The in.vent Panels of Normal Pregnancy – consecutive donations are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments, mainly for study schemes regarding cut offs and studies for certification. The panel consists of 30 samples from 10 different donors. 1 sample per donor per trimester is provided.

Product No.:	FPPA 4.7.2.1
Number of Donors	10
Consecutive Donations	3 (1 per trimester)
Number of Samples	30
Volume per Samples	1.0 ml, different volumes on request
Source Material	true human serum off the clot
Details given	product number, sample number, matrix, volume, age, gender, week of pregnancy, minimum shelf life, scheme of processing, results for several laboratory parameters on availability
Materials Safety	All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	680.00 €/Panel

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
inquiry.inventdiagnostica.de



made in germany

www.inventdiagnostica.de