



Product Datasheet

In.QualiSet SARS-CoV-2 Validation Panel

True Human Serum off the Clot

Application:

The In.QualiSet SARS-CoV-2 Validation Panel is designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. This In.QualiSet Panel contains SARS-CoV-2 positive serum samples and negative serum samples with 16 IgG titer of diseases, which can interfere with SARS-Cov-2 IgG, IgM and IgA antibody tests.

Product No.:	DSPA 4.1.9.8.1
Number of Samples	96
Volume per Sample	150 µl
Matrix	true human serum off the clot, without additives
Details given	product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life, date of blood donation
Content	<ul style="list-style-type: none"> • 6 x SARS-CoV-2 antibody positive samples <ul style="list-style-type: none"> ○ 4 donors (2 donors with consecutive samples) with previous positive SARS-CoV-2 PCR test • 90 x human samples (pre-SARS-CoV-2) measured positive on: <ul style="list-style-type: none"> ○ 81 x Adenovirus IgG ○ 66 x Cytomegalie Virus (CMV) IgG ○ 73 x Herpes Simplex 1+2 IgG ○ 86 x Influenzaviruses Type A IgG ○ 72 x Influenzaviruses Type B IgG ○ 73 x Mumps Virus (Parotitis) IgG ○ 20 x Mycoplasma pneumoniae IgG ○ 90 x Parainfluenza 1/2/3 IgG ○ 50 x Parvovirus B19 IgG ○ 89 x Respiratory Syncytial Virus IgG ○ 87 x Varicella Zoster Virus IgG ○ 87 x Anti-Epstein-Barr-Virus-Capsid-Antigen IgG ○ 18 x Anti-Enteroviruses IgG ○ 56 x Anti-Measles Virus IgG ○ 1 x Anti-Hanta-Virus Pool 1 "Eurasia" IgG ○ 45 x Haemophilus influenzae B IgG
Infectious Status	All pre-SARS-CoV-2 samples are tested PCR negative for HIV, HBV, HCV.
Intended Use	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -18°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de