



## Product Datasheet

### Quintet of Donors >25 years

#### Serum and plasma

#### Application:

Quintets of Preanalytics are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

<b>Product No.:</b>	NDSA 4.3.1.4
<b>Number of Samples</b>	5
<b>Volume per Sample</b>	1.0 ml
<b>Source Material</b>	serum and plasma of self-declared healthy donors
<b>Details given</b>	product number, volume, matrix, minimum shelf life, medication and contraceptives, gender, age, assay system, smoker/non-smoker
<b>Donors</b>	<ul style="list-style-type: none"> <li>• age: &gt; 25 years</li> <li>• no diabetes</li> <li>• BMI &lt;30</li> </ul>
<b>Content of Quintet</b>	<ul style="list-style-type: none"> <li>• 2 x 1.0 ml serum</li> <li>• 1 x 1.0 ml K2 EDTA plasma</li> <li>• 1 x 1.0 ml K3 EDTA plasma</li> <li>• 1 x 1.0 ml Li-Heparin plasma</li> <li>• 1 x 1.0 ml Na-Citrat plasma</li> </ul>
<b>Disease State Testing</b>	on request
<b>Materials Safety</b>	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	24 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](http://inquiry.inventdiagnostica.de)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)