



Product Datasheet

Quintet of Normal Donors

Serum and plasma

Application:

Quintets of Preanalytics are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

Product No.:	NDSA 4.3.1.2
Number of Samples	5
Volume per Sample	1.0 ml
Source Material	serum and plasma of self-declared healthy donors
Details given	product number, volume, matrix, minimum shelf life, medication and contraceptives, gender, age, assay system, smoker/non-smoker
Donors	<ul style="list-style-type: none"> • age: 25-65 years • non-smokers • no diabetes • BMI <30
Content of Quintet	<ul style="list-style-type: none"> • 2 x 1.0 ml serum • 1 x 1.0 ml K2 EDTA plasma • 1 x 1.0 ml K3 EDTA plasma • 1 x 1.0 ml Li-Heparin plasma • 1 x 1.0 ml Na-Citrat plasma
Disease State Testing	on request
Materials Safety	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de