

Product Datasheet

SARS-CoV-2 Distinction Panel

True Human Serum off the Clot

Application:

The in.vent SARS-CoV-2 Distinction Panel is designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. The entire Range of Laboratory Test Results is represented within the Panel. The Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

Product No.:	DSPA 4.1.9.1.2
Number of Samples	90
Volume per Sample	0.2-1.5 ml
Source Material	true human serum off the clot
Details given	product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life, date of blood donation
Diagnoses	10 samples of Chl. Pneumoniae 10 samples of CMV 10 samples of EBV 10 samples of HCV 10 samples of HSV 10 samples of Mycopl. Pneumoniae 10 samples of Rheumatoid Factor 10 samples of Borrelia 10 samples of ANA
Materials Safety	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -18°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	6,000.00 €/Panel

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de