



Product Datasheet

anti-TSH receptor (TRAb) Pool

Pool of Human Plasma, TRAb positive

Application:

Pool of Human Plasma, TRAb positive, is used for production purposes like production of standards, calibrators or controls for IVDs or EQAS.

Product No.:	DSPO 4.1.6.4.5
Volume	5.0 l, individual volumes on request
Source Material	native human sodium-citrate plasma, pooled
Titer	TRAb \geq 1000 IU/ml (Siemens Immulite)
Details given	titer (TRAb, anti-TPO, anti-Tg), ratio (TRAb/anti-Tg/anti-TPO), volume, matrix, test kit, test manufacturer
Ratio anti-Tg/anti-TPO, TRAb	on availability
Filtration	0.2 μ m
Bottling	1.0 l
Test Samples	available
Disease State Testing	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV RNA, Syphilis.
Materials Safety	<p>These methods cannot eliminate possible infectivity of the material. All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material.</p> <p>All products are for research or manufacturing purposes only and are not for use in human therapeutic procedures.</p>
Shipping	on dry ice
Storage	\leq -15°C
Shelf Life	48 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
inquiry.inventdiagnostica.de



made in germany

www.inventdiagnostica.de