



## Product Datasheet

### Rheumatoid Arthritis Control

### True Human Serum off the Clot

#### Application:

in.vent Quality Controls are designed to be used in Laboratory Quality Control Settings, in Quality Departments and R&D Units. They are applicable for solving tasks of precision, accuracy or correctness of in-vitro diagnostic assay measurements and analyte components throughout the whole value chain of in-vitro diagnostic test kits or assays.

<b>Product No.:</b>	QC 2.1.2
<b>Source Material</b>	true human serum off the clot
<b>Bottling</b>	1.0 ml each vial
<b>Parameters</b>	Rheumatoid Factor IgG, Rheumatoid Factor IgM, Rheumatoid Factor IgA, anti-CCP
<b>Provided Levels</b>	2 levels positive concentrations of RF IgG + IgM + IgA and anti-CCP
<b>Stability</b>	5 years at $\leq 15$ °C, 1 year at 2-8 °C
<b>Form</b>	liquid
<b>Disease State Testing</b>	tested by CE marked test systems and found negative for HIV I/II antibodies, HBs antigen, HCV antibodies
<b>Pre-Samples</b>	available
<b>Delivered Information on Certificate</b>	product number, matrix, processing, homogeneity, target values, reconstitution instruction
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)