



Product Datasheet

Rheumatoid Arthritis Pool

True Human Serum off the Clot

Application:

Pool of True Human Serum off the Clot, Rheumatoid Arthritis, is used for production purposes like production of standards, calibrators or controls for IVDs or EQAS.

Product No.:	DSPO 4.1.1.1.1
Volume	1.0 l individual volumes on request
Source Material	true human serum off the clot
Biomarker	Rheumatoid Factor
Details given	titer, volume, matrix, test kit, test manufacturer
Bottling	100.0 ml
Test Samples	available
Disease State Testing	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV RNA, HBV DNA, HCV, RNA, Syphilis.
Materials Safety	All material of human origin is to be regarded as potentially infectious. Therefore standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	48 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de