



Product Datasheet

Parainfluenza 1/2/3 IgG positive Panel

True Human Serum off the Clot

Application:

The in.vent Panels of Parainfluenza 1/2/3 IgG positive are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

Product No.:	DSPA 4.1.4.1.24
Number of Samples	10
Volume	1.0 ml each sample
Matrix	true human serum off the clot
Details given	product number, sample number, volume, matrix, age, gender, measured laboratory result, assay system, minimum shelf life
Parameter	Parainfluenza 1/2/3 IgG (titers, further parameters and diagnoses on availability)
Disease State Testing	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV, HIV I RNA, HCV, RNA, HBV DNA, HCV RNA, Syphilis.
Materials Safety	All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de