



Product Datasheet

Thyroglobulin (Tg) free Panel

True Human Serum off the Clot

Application:

The in.vent Tg free Panels are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

| | |
|-------------------------------------|--|
| Product No.: | DSPA 4.1.6.7 |
| Number of Samples | 10 |
| Volume | 1.0 ml each sample |
| Matrix | true human serum off the clot |
| Details given | product number, sample number, volume, matrix, age, gender, measured laboratory result, assay system, minimum shelf life |
| Parameter | Tg free (titers and further parameters on availability): <ul style="list-style-type: none">• 2 samples measured negative• 8 samples measured positive |
| Materials Safety | All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals. |
| Storage | ≤ -15°C |
| Minimum Shelf Life | 24 months |
| Manufacturing Quality System | in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485. |
| Pricing | Please inquire for your individual quote. |

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
[inquiry.inventdiagnostica.de](https://www.inventdiagnostica.de/inquiry)



made in germany

www.inventdiagnostica.de