



## Product Datasheet

### Diabetes Panel with selected Medication

#### True Human Serum off the Clot

#### Application:

The in.vent Diabetes Panels with selected Medication are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments.

<b>Product No.:</b>	DSPA 4.1.3.2.3
<b>Number of Samples</b>	20
<b>Volume</b>	1.0 ml each sample
<b>Matrix</b>	true human serum off the clot, optional: plasma
<b>Details given</b>	product number, sample number, volume, matrix, age, gender, measured laboratory result, assay system, minimum shelf life
<b>Diagnoses</b>	<ul style="list-style-type: none"> <li>• 10 samples with diagnosis diabetes type I with typical medication, optional: correlating parameters</li> <li>• 10 samples with diagnosis diabetes type II with typical medication</li> </ul>
<b>Materials Safety</b>	<p>All material of human origin is to be regarded as potentially infectious. Therefore, standard precautions are necessary, if handling the material.</p> <p>All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.</p>
<b>Storage</b>	≤ -15°C
<b>Minimum Shelf Life</b>	24 months
<b>Manufacturing Quality System</b>	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](http://inquiry.inventdiagnostica.de)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)