



Product Datasheet

In.QualiSet Interference Panel

True Human Serum off the Clot, Plasma

Application:

The In.QualiSet Interference Panels are designed for the use throughout the whole value chain of in-vitro diagnostic test kits, e.g. for R&D and Quality Control Departments. The entire Range of Laboratory Test Results is represented within the Panels. The Samples consist of native Human Bio-Material with no Additives. Titers are based upon Availabilities.

Product No.:	IFPA 4.5.1.1
Number of Samples	96
Volume per Sample	150 µl
Source Material	true human serum off the clot, K2 EDTA plasma, K3 EDTA plasma, Li-Heparin plasma, Na-Citrate plasma
Details given	product number, sample number, matrix, volume, age, gender, measured laboratory result, assay system, minimum shelf life
Diagnoses	<ul style="list-style-type: none"> • 8 x HAMA positive serum • 6 x Rheumatoid Factor Antibodies • 2 x anti-E. coli Antibodies • 3 x Vitamin D3 • 6 x K3 EDTA plasma • 3 x K2 EDTA plasma • 6 x Li-Heparin plasma • 6 x Na-Citrate plasma • 5 x Ibuprofen • 5 x Aspirin • 4 x Diclofenac • 5 x Insulin • 5 x Antihypertensives • 6 x Bilirubin (icteric) • 2 x Haemoglobin (haemolytic) • 8 x Triglyceride (lipaemic) • 8 x Progesterone • 8 x Testosterone
Materials Safety	All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
Storage	≤ -15°C
Minimum Shelf Life	24 months
Manufacturing Quality System	in.vent Diagnostica GmbH operates a quality management system, which complies with the requirements of DIN EN ISO 9001 and DIN EN ISO 13485.
Pricing	Please inquire for your individual quote.

Contact & Ordering Information



Contact for your Order
+49 (0)3302 55 199 250
sales@inventdiagnostica.de

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

Send us your individual Request:
inquiry.inventdiagnostica.de



made in germany

www.inventdiagnostica.de