



## Product Datasheet

### invent Study Site: Whole Blood Donations

#### Whole Blood (K3EDTA)

#### Application:

invent runs a company-owned and fully equipped study site. Together with our well-trained staff we conduct on demand donation days and are in full control over all steps of pre-analytics.

<b>Product No.:</b>	SV 7.1.1.2
<b>Source Material</b>	normal donors, self-declared healthy, fresh frozen
<b>Details given</b>	product number, matrix, preanalytics, processing, minimum shelf life
<b>Matrix</b>	whole blood with K3EDTA
<b>Volume</b>	~500ml per donation
<b>Pre-analytic</b>	donor specifications: gender, age processing: frozen within 24h
<b>Materials safety</b>	These methods cannot eliminate possible infectivity of the material. All material of human origin is to be regarded as potentially infectious. Therefore standard precautions are necessary, if handling the material. All products are for research or manufacturing purposes only and are not for direct therapeutic or diagnostic use in humans or animals.
<b>Storage</b>	< -15 °C
<b>Delivery</b>	fresh frozen on dry ice, delivered in between 72 hours after blood collection
<b>Infectious Disease Testing</b>	Each individual donation has been tested with methods cleared in compliance with European Directive 98/79/EC and found negative for HBsAg, anti-HIV I+II, anti-HCV
<b>Pricing</b>	Please inquire for your individual quote.

# Contact & Ordering Information



**Contact for your Order**  
**+49 (0)3302 55 199 250**  
**[sales@inventdiagnostica.de](mailto:sales@inventdiagnostica.de)**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

**Send us your individual Request:**  
**[inquiry.inventdiagnostica.de](http://inquiry.inventdiagnostica.de)**



made in germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)