

*We are human*

**in.vent**  
D I A G N O S T I C A G M B H

## External Quality Assessment Scheme and CE-Controls



**in.vent** Diagnostica GmbH  
Neuendorfstraße 17  
16761 Hennigsdorf  
Germany

[www.inventdiagnostica.de](http://www.inventdiagnostica.de)



EQAS

For more than **20 years**, **in.vent** has been developing and producing interdisciplinary quality controls specifically tailored to the requirements of our partners. We work closely with international proficiency test providers and also manufacture OEM controls. Development of customized CE-labelled controls is also available. All quality controls are based on preferably native human biomaterials from ethically approved collections. Our focus is on the highest quality, reproducibility and practical relevance.

Thanks to our many years of experience, we offer reliable solutions for a wide range of diagnostic areas. The long-standing loyalty of our customers is an expression of their trust in our quality standards and underlines their high level of satisfaction with our products and services.

1. Alcohol Abuse

| Type of Control                                  | Product No. | Analyte                               | Method    | Sample                             | Minimum Shelf Life |
|--------------------------------------------------|-------------|---------------------------------------|-----------|------------------------------------|--------------------|
| Carbohydrate deficient transferrin (CDT) Control | 840000510   | Asialotransferrin, disialotransferrin | HPLC      | 1.0mL serum, lyophilized           | 2 years            |
| Ethyl Glucuronide (EtG) Control                  | 840000511   | Ethyl glucuronide                     | LC-MS/MS* | 1.0mL urine, lyophilized           | 2 years            |
| NEW Phosphatidylethanol (PEth) Control           | 840000512   | Phosphatidylethanol                   | HPLC      | 0.5-1.0mL whole blood, lyophilized | 1 year             |

2. Allergy

| Type of Control            | Product No. | Analyte                                                              | Method           | Sample                             | Minimum Shelf Life |
|----------------------------|-------------|----------------------------------------------------------------------|------------------|------------------------------------|--------------------|
| Inhalation Allergy Control | 840000500   | e.g. d01, d02, e01, g06, g12, t03, t04, w06, total IgE (*)           | FEIA** ImmunoCAP | 1.0mL serum or plasma, lyophilized | 2 years            |
| Insect Allergy Control     | 840000501   | e.g. i01, i03, total IgE (*)                                         | FEIA** ImmunoCAP | 1.0mL serum or plasma, lyophilized | 2 years            |
| Nutrition Allergy Control  | 840000502   | e.g. f02, f03, f04, f05, f09, f13, f14, f17, f49, f84, total IgE (*) | FEIA** ImmunoCAP | 1.0mL serum or plasma, lyophilized | 2 years            |

\*Different allergen compositions available  
\*\*approved for all common test systems

3. Autoimmunity

| Type of Control                                             | Product No. | Analyte                                                                                                                                                                                                                                                                                                          | Method                     | Sample                             | Minimum Shelf Life             |
|-------------------------------------------------------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------|------------------------------------|--------------------------------|
| Anti-Extractable nuclear antigens (ENA) Control*            | 840000400   | Anti-Jo 1 (histidyl-tRNA synthetase), anti-PM/Scl, anti-SS-A/Ro, anti-SS-A/Ro (52 kD proteins), anti-SS-A/Ro (60 kD - proteins), anti-SS-B/La, anti-Sm, anti-SmD, anti-topoisomerase-1 (Scl 70), anti-U1 snRNP                                                                                                   | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Antinuclear antibodies (ANA) Control*                       | 840000402   | Autoantibodies against cell nuclei (ANA), autoantibodies against dsDNA, collagenosis screening, anti-DFS70/LEDGF                                                                                                                                                                                                 | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Antiphospholipid syndrome (APS) Control*                    | 840000403   | Anti-cardiolipin (IgA), anti-cardiolipin (IgG), anti-cardiolipin (IgM), anti-β2-glycoprotein (IgA), anti-β2-glycoprotein (IgG), anti-β2-glycoprotein (IgM)                                                                                                                                                       | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Autoimmune dermatosis Control*                              | 840000404   | Anti-BPAG 180, anti-BPAG 230, anti-desmoglein 1, anti-desmoglein 3, anti-desmoplakin I/II, anti-spiny cell desmosomes, anti-epidermal basement membrane                                                                                                                                                          | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Autoimmune endocrinopathy Control*                          | 840000405   | Anti-adrenocortical, anti-P450 C17-hydroxylase, anti-P450 C21-hydroxylase, anti-P450 ssc-hydroxylase                                                                                                                                                                                                             | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Diabetes mellitus Control*                                  | 840000406   | Anti-glutamate decarboxylase (GAD), anti-islet cell (ICA), anti-insulin (IAA), anti-tyrosine phosphatase (IA2), anti-zinc transporter (ZnT8A)                                                                                                                                                                    | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Gluten sensitive enteropathy Control                        | 840000407   | Anti-endomysium IgA, anti-gliadin IgA, anti-gliadin IgG, anti-transglutaminase IgA, anti-deamidated gliadin IgA, anti-deamidated gliadin IgG                                                                                                                                                                     | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Hepatopathy Control*                                        | 840000408   | Anti-actin, anti-H+/K+-ATPase, anti-microsomes (LKM 1), anti-mitochondria, anti-parietal cells, anti-smooth muscles                                                                                                                                                                                              | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Myasthenia gravis Control*                                  | 840000409   | Anti-acetylcholine receptors, anti-muscle-specific receptor tyrosine kinase (MuSK), anti-skeletal muscle, anti-titin                                                                                                                                                                                             | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Paraneoplastic Neuropathy Control*                          | 840000410   | Anti-glial cell nuclei antibodies (AGNA), anti-amphiphysin 1, anti-HuD, anti-neuron nuclei, anti-Purkinje cells, anti-Ri/Nova 1, anti-SOX1, anti-Ta (Ma 2), anti-Yo (CDR62)                                                                                                                                      | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Peripheral neuropathy Control*                              | 840000411   | anti-asialoganglioside (AGA1), anti-disialoganglioside (GD1a), anti-disialoganglioside (GD1b), anti-monosialoganglioside (GM1), anti-monosialoganglioside (GM2), anti-monosialoganglioside (GM3), anti-myelin-associated glycoprotein (MAG), anti-quadrosialoganglioside (GQ1b), anti-trisialoganglioside (GT1b) | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Rheumatoid Arthritis Control*                               | 840000412   | Rheumatoid factor, rheumatoid factor isotypes, anti-CCP (cyclic citrulline peptide), anti-MCV (mutant)                                                                                                                                                                                                           | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| Vasculitis and Glomerulonephritis Control*                  | 840000413   | Anti-MPO, anti-PR3, anti-glomerular basement membrane (GBM), c-ANCA, p-ANCA                                                                                                                                                                                                                                      | CIA, IFT, Blot, ELISA, RIA | 0.3mL serum/plasma, liquid, frozen | 2 years                        |
| NEW Anti-Saccharomyces cerevisiae Antibodies (ASCA) Control | 840000580   | ASCA IgG<br>ASCA IgA                                                                                                                                                                                                                                                                                             | ELISA                      | 0.5 ml serum lyophilised           | upon request from the customer |

\*approved for all common test systems

4. Endocrinology

| Type of Control                              | Product No. | Analyte                                                                                                                                                                                         | Method                                                    | Sample                    | Minimum Shelf Life             |
|----------------------------------------------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------|---------------------------|--------------------------------|
| Adrenal & Pituitary gland Control*           | 840000440   | 1-39 ACTH, DHEA, DHEA-S, HGH, prolactin                                                                                                                                                         | ECLIA, RIA, CLEIA                                         | 3mL serum, lyophilized    | 2 years                        |
| Aldosterone Control                          | 840000441   | Aldosterone                                                                                                                                                                                     | CLIA                                                      | 1mL serum, lyophilized    | 2 years                        |
| Bone Metabolism Control*                     | 840000442   | Osteocalcin, CTx (CrossLaps) Bone-specific alkaline phosphatase (BAP), Procollagen type 1 amino-terminal propeptide (P1NP)                                                                      | Kryptor, Liaison XL, Cobas Elecsys, Modular E170* (ECLIA) | 2mL serum, lyophilized    | 2 years                        |
| Calcium Metabolism & Special thyroid Control | 840000443   | Parathyroid hormone 1-84 (PTH 1-84), calcitonin (hCT), thyroglobulin (Tg)                                                                                                                       | Kryptor, CLIA, ECLIA*                                     | 3mL serum, lyophilized    | 2 years                        |
| Cortisol Salivary Control                    | 840000444   | Cortisol                                                                                                                                                                                        | ECLIA, ELISA                                              | 1mL serum, lyophilized    | 2 years                        |
| Fertility Control*                           | 840000445   | Anti-Müller-Hormone (AMH), Follicle-stimulating hormone (FSH), Luteinizing hormone (LH)                                                                                                         | ECLIA                                                     | 3mL serum, lyophilized    | 2 years                        |
| Hematopoietic Growth Factors Control         | 840000446   | Erythropoietin (EPO), Thrombopoietin (THPO)                                                                                                                                                     | ELISA                                                     | 1mL serum, lyophilized    | 2 years                        |
| Pancreatic & Gastrointestinal Control*       | 840000447   | C-peptide, Gastrin, Glucagon, Insulin, Proinsulin                                                                                                                                               | ELISA, ECLIA, CLIA, RIA                                   | 3mL serum, lyophilized    | 2 years                        |
| Thyroid Autoantibodies Control*              | 840000448   | TSH receptor antibody (TRAK), Thyroid peroxidase antibody (MAK), Thyroglobulin antibody (TAK)                                                                                                   | ECLIA                                                     | 1mL serum, lyophilized    | 2 years                        |
| Thyroid Hormones Control*                    | 840000449   | Triiodothyronine (Total-T3), thyroxine (total-T4), TSH, f-triiodothyronine (free-T3), f-thyroxine (free-T4)                                                                                     | ECLIA                                                     | 3mL serum, lyophilized    | 2 years                        |
| Virilization Control                         | 840000450   | 17-Hydroxy Progesterone (17 OH-P), Androstenedione (ASD)                                                                                                                                        | ECLIA, CLIA                                               | 1mL serum, lyophilized    | 2 years                        |
| NEW Bile acid control                        | 840000573   | cholic acid, chenodeoxycholic acid                                                                                                                                                              | LC-MS/MS                                                  | 1.0 mL serum, lyophilised | upon request from the customer |
| NEW Thyroid Multicontrol 1-physiological     | 840000575   | Anti-TPO < 25.0 IU/ml<br>Anti Tg < 25.0 IU/ml<br>TRAK 0.5 -1.2 IU/ml<br>fT3 3.0 -6.0 pmol/L<br>fT4 13.0 - 16.0 pmol/L<br>T3 1.0 - 3.0 nmol/L<br>T4 70.0 - 140.0 nmol/L<br>TSH 1.0 - 3.0 mU/L    | ECLIA, CMIA                                               | 3.0 mL serum, lyophilised | upon request from the customer |
| NEW Thyroid Multicontrol 2-Hyperthyroidism   | 840000576   | Anti-TPO < 25.0 IU/ml<br>Anti Tg <25.0 IU/ml<br>TRAK 5.0 -8.0 IU/ml<br>fT3 12.0 - 18.0 pmol/L<br>fT4 25.0 - 32.0 pmol/L<br>T3 3.0 - 8.0 nmol/L<br>T4 170.0 - 250.0 nmol/L<br>TSH 0.1 - 0.8 mU/L | ECLIA, CMIA                                               | 3.0 mL serum, lyophilised | upon request from the customer |
| NEW Thyroid Multicontrol 3-Hypothyroidism    | 840000577   | Anti-TPO > 250.0 IU/ml<br>Anti Tg > 250.0 IU/ml<br>TRAK 0.5 -1.2 IU/ml<br>fT3 1.0 -3.0 pmol/L<br>fT4 5.0 - 12.0 pmol/L<br>T3 0.2 - 1.2 nmol/L<br>T4 30.0 - 80.0 nmol/L<br>TSH 5.0 - 10.0 mU/L   | ECLIA, CMIA                                               | 3.0 mL serum, lyophilised | upon request from the customer |
| NEW Glucagon Control                         | 840000578   | Glucagon                                                                                                                                                                                        | ELISA                                                     | 1.0 mL serum, lyophilised | upon request from the customer |
| NEW Calcitonin Control                       | 840000579   | Calcitonin                                                                                                                                                                                      | CLIA                                                      | 1.0 mL serum, lyophilised | upon request from the customer |

\*approved for all common test systems

5. Gastroenterology

| Type of Control                                       | Product No. | Analyte                                                                                                            | Method                                | Sample                   | Minimum Shelf Life |
|-------------------------------------------------------|-------------|--------------------------------------------------------------------------------------------------------------------|---------------------------------------|--------------------------|--------------------|
| Calprotectin & Elastase faecal Control                | 840000470   | Calprotectin, Calprotectin (qualitative), pancreatic Elastase, pancreatic Elastase (qualitative)                   | CLIA, EIA, ELISA, Lateral flow assay* | 0.5g faeces, lyophilized | 1 year             |
| Helicobacter pylori antigen faecal Control            | 840000471   | Helicobacter pylori antigen                                                                                        | CLIA, Rapid Test                      | 0.5g faeces, lyophilized | 2 years            |
| Hemoglobin (Hb) - Haptoglobine complex faecal Control | 840000472   | Hemoglobin, hemoglobin (qualitative), hemoglobin-haptoglobin complex, hemoglobin-haptoglobin complex (qualitative) | ELISA, CLIA                           | 0.5g faeces, lyophilized | 2 years            |
| Occult blood faecal Control                           | 840000473   | Hemoglobin in the stool                                                                                            | Turbidimetry, IFOBT                   | 0.5g faeces, lyophilized | 6 months           |

\*approved for all common test systems

6. Urology

| Type of Control                              | Product No. | Analyte                                                                                                                                                                                                                                                     | Method                                                                                                  | Sample                     | Minimum Shelf Life             |
|----------------------------------------------|-------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------|----------------------------|--------------------------------|
| Urine Standard profile, qualitative Control  | 840000530   | Urine albumin, ascorbic acid, bilirubin, protein, erythrocytes, glucose, ketone bodies, leukocytes, nitrite, urobilinogen, hCG, pH value, specific gravity, Crystals                                                                                        | Clintek (Siemens) turbidimetry various strip test immunoassay                                           | 10mL urine, liquid, frozen | 2 years                        |
| Urine Standard profile, quantitative Control | 840000531   | Protein, Glucose, Uric acid, Urea, Creatinine, Microalbumin, Calcium, Chloride Iron, Potassium, Copper, Magnesium,Sodium, Phosphate, alpha-amylase, Osmolality, Pancreas Amylase, alpha-1 microglobulin, beta 2 microglobulin, IgG, Transferrin, hCG intact | Roche Cobas 8000, Agilent Cary 60, AAS NovAA, BN ProSpec Siemens, Siemens BN2, Abbott Architect (ECLIA) | 3mL urine, lyophilized     | 1 year                         |
| NEW Cotenin in Urine Control                 | 840000570   | Cotenin                                                                                                                                                                                                                                                     | Immunoassay / LC-MS/MS                                                                                  | 1.0 ml urine, lyophilised  | upon request from the customer |
| NEW Free Cortisol in Urine Control           | 840000571   | Cortisol                                                                                                                                                                                                                                                    | ECLIA                                                                                                   | 1.0 ml urine, lyophilised  | upon request from the customer |

7. Immunology

| Type of Control                            | Product No. | Analyte                                                                                                                                              | Method                                        | Sample                                 | Minimum Shelf Life             |
|--------------------------------------------|-------------|------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|----------------------------------------|--------------------------------|
| Gammopathies Control*                      | 840000480   | Free light chains type kappa (quantitative), free light chains type lambda (quantitative), gammopathie types, IgA, IgG, IgM, kappa/lambda (quotient) | Immunonephelometry, Capillary electrophoresis | 1mL serum, lyophilized                 | 2 years                        |
| IgG subclasses Control*                    | 840000481   | IgG1, IgG2, IgG3, IgG4                                                                                                                               | Immuno-nephelometry                           | 1mL serum, lyophilized                 | 2 years                        |
| Immunophenotyping Control                  | 840000482   | Analytes according to case description, T, NK and B cell subpopulations including naïve, memory, effector, activated, TH1/2 and regulatory cells     | flow cytometry                                | EDTA whole blood not stabilized        | 24 hours                       |
| Interferon Gamma Release Control*          | 840000483   | TB1 NIL, TB2 NIL, MITO IFN                                                                                                                           | Quantiferon TB Gold Plus                      | Li-heparin whole blood, not stabilized | 24 hours                       |
| NEW Basophil Activation Test Control (BAT) | 840000574   | Basophil activation by honeybee and wasp (yellow jacket) allergens                                                                                   | Flow Cytometer                                | 4.5 ml EDTA whole blood                | upon request from the customer |

\*approved for all common test systems

8. Pregnancy Diagnostics

| Type of Control                  | Product No. | Analyte                                                                     | Method     | Sample                        | Minimum Shelf Life |
|----------------------------------|-------------|-----------------------------------------------------------------------------|------------|-------------------------------|--------------------|
| 1st trimester pregnancy Control* | 840000490   | Free beta-hCG, PAPP-A                                                       | FIA, ECLIA | 1mL serum, lyophilized        | 1 year             |
| 2nd trimester pregnancy Control* | 840000491   | AFP, hCG intact, hCG intact + beta hCG                                      | FIA, ECLIA | 1mL serum, lyophilized        | 1 year             |
| Preeclampsia Control*            | 840000492   | Placental growth factor (PLGF), Soluble fms-like tyrosine kinase-1 (sFLT-1) | FIA, ECLIA | 1mL serum, lyophilized        | 2 years            |
| Fetal RhD-Factor Control         | 840000590   | Fetal RhD factor                                                            | qPCR       | 1 - 2mL plasma liquid, frozen | 1 year             |

\*approved for all common test systems

9. Medical Drugs

| Type of Control                                             | Product No. | Analyte                                                                                                         | Method                                          | Sample                              | Minimum Shelf Life |
|-------------------------------------------------------------|-------------|-----------------------------------------------------------------------------------------------------------------|-------------------------------------------------|-------------------------------------|--------------------|
| Antibiotics, Cytostatics & Cardiac glycosides Control*      | 840000460   | Amikacin, Digitoxin, Gentamicin, Methotrexate, Tobramycin, Vancomycin                                           | ECLIA, PETINIA, Fluorescence polarization, CMIA | 5mL serum, lyophilized              | 2 years            |
| Antiepileptics, Cardiac Glycosides & Antiasthmatic Control* | 840000462   | Carbamazepine, Digoxin, Ethosuximide, Lithium, Phenobarbital, Phenytoin, Primidone, Theophylline, Valproic acid | HPLC, PETINIA, EIA, ECLIA, CHROM, KIMS          | 5mL serum, lyophilized              | 2 years            |
| Paracetamol Control*                                        | 840000464   | Paracetamol                                                                                                     | CEDIA (Roche)                                   | 1mL serum, lyophilized              | 2 years            |
| Salicylic acid Control*                                     | 840000465   | Salicylic acid                                                                                                  | FPIA (Roche)                                    | 1mL serum, lyophilized              | 2 years            |
| Immunosuppressants Control                                  | 840000553   | Cyclosporine, Sirolimus, Everolimus, Tacrolimus                                                                 | LC-MS/MS                                        | 2 mL human whole Blood, lyophilised | 1 year             |

\*approved for all common test systems

10. Liquor Diagnostics

| Type of Control                                    | Product No. | Analyte                                                                                                                                                                                                             | Method | Sample                         | Minimum Shelf Life |
|----------------------------------------------------|-------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------|--------------------------------|--------------------|
| Cerebrospinal Fluid (CSF) Standard Profile Control | 840000560   | Glucose quotient (cerebrospinal fluid/serum), glucose quotient - diagnostic interpretation, glucose in the cerebrospinal fluid, interpretation of the findings report, lactate, lactate - diagnostic interpretation | ECLIA  | 0.5mL CSF + serum, lyophilised | 9 months           |

11. Infection Diagnostics

| Type of Control             | Product No. | Analyte            | Method | Sample                    | Minimum Shelf Life             |
|-----------------------------|-------------|--------------------|--------|---------------------------|--------------------------------|
| Procalcitonin (PCT) Control | 840000520   | PCT                | ECLIA  | 1mL serum, lyophilized    | 2 years                        |
| Anti TBE Control            | 840000572   | TBE IgG<br>TBE IgM | ELISA  | 0.5 ml serum, lyophilised | upon request from the customer |

12. Oncology

| Type of Control                           | Product No. | Analyte  | Method | Sample                 | Minimum Shelf Life |
|-------------------------------------------|-------------|----------|--------|------------------------|--------------------|
| Prostate-specific -antigen (PSA) Control* | 840000540   | PSA/fPSA | ECLIA  | 1mL serum, lyophilized | 2 years            |

\*approved for all common test systems

13. CE Controls

| Type of Control                                                    | Product No. | Analyte                                                                                                                   | Method                                                                                                                  | Sample                         | Minimum Shelf Life |
|--------------------------------------------------------------------|-------------|---------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|--------------------------------|--------------------|
| In.QualiCheck Calprotectin Control 1*<br>QC 2.1.20                 | 810000300   | Faecal calprotectin < 50µg/g                                                                                              | Bühlmann Calprotectin fCal turbo Kit, (PETIA) BA 200 (Bio-Systems)                                                      | 0.5g Human faeces, lyophilized | 2 years            |
| In.QualiCheck Calprotectin Control 2*<br>QC 2.1.21                 | 810000301   | Faecal calprotectin 50 - 200µg/g                                                                                          | Bühlmann Calprotectin fCal turbo Kit, (PETIA) BA 200 (Bio-Systems)                                                      | 0.5g Human faeces, lyophilized | 2 years            |
| In.QualiCheck Calprotectin Control 3*<br>QC 2.1.22                 | 810000302   | Faecal calprotectin > 200µg/g                                                                                             | Bühlmann Calprotectin fCal turbo Kit, (PETIA) BA 200 (Bio-Systems)                                                      | 0.5g Human faeces, lyophilized | 2 years            |
| In.QualiCheck Prostate Specific Antigen (PSA) Controls<br>QC 2.1.8 | 810000306   | Level 1<br>PSA, total 2 - 6ng/mL<br>PSA, free 1 - 3ng/mL<br><br>Level 2<br>PSA, total 6 - 10ng/mL<br>PSA, free 3 - 6ng/mL | PSA Cancer Check, POCT, concile® Ω100 or Roche Elecsys®, ECLIA, Cobas 8000 e602 or Siemens Atellica®, CLIA, IM Analyzer | 0.5mL Human serum, lyophilized | 2 years            |

\*Available through April 2027

# Contact & Ordering Information



## General Contact

**+49 (0)3302 55 199 0**

**info@inventdiagnostica.de**

Neuendorfstraße 17 | 16761 Hennigsdorf | Germany

---

## Contact for your Order

**sales@inventdiagnostica.de**

**Please send us your individual requests via:  
inquiry.inventdiagnostica.de**



Version September 2026



Made in Germany

**www.inventdiagnostica.de**